

# GeneFinder™ COVID-19 Fast RealAmp Kit

**REF** IFMR-49 **Σ** 100 tests/Kit

**Store at -20°C or below.**  
**Shelf life is 20 months after manufacturing**

**CE-IVD**  
**For ABI7500 Fast, CFX96, QuantStudio 5,  
qTOWER<sup>3</sup> and Gentier 96E**

## INTENDED USE

GeneFinder™ COVID-19 Fast RealAmp Kit is designed for qualitative detection of SARS-CoV-2 RNA with real time reverse transcription PCR from Nasopharyngeal & Oropharyngeal (Throat) swab, Sputum and Bronchoalveolar Lavage Fluid (BAL) samples within 50 minutes.

## KIT COMPONENT

COVID-19 Fast RealAmp Kit	100 tests / Kit
COVID-19 Fast Reaction Buffer <sup>1)</sup>	400 µL
COVID-19 Fast Enzyme Mixture <sup>2)</sup>	200 µL
COVID-19 Fast Probe Mixture <sup>3)</sup>	500 µL
COVID-19 Fast Positive Control <sup>4)</sup>	50 µL
COVID-19 Fast Negative Control <sup>5)</sup>	50 µL
DEPC Water <sup>6)</sup>	400 µL

<sup>1)</sup> COVID-19 Fast Reaction Buffer: Containing Tris-Cl, MgCl<sub>2</sub> and dNTPs etc.

<sup>2)</sup> COVID-19 Fast Enzyme Mixture: Containing RTase, Taq polymerase, and RNase inhibitor etc.

<sup>3)</sup> COVID-19 Fast Probe Mixture: Primer pairs and Probes for amplification and detection of each target.

<sup>4)</sup> COVID-19 Fast Positive Control; Positive control for targets.

<sup>5)</sup> COVID-19 Fast Negative Control; Ultrapure quality water, PCR-grade.

<sup>6)</sup> DEPC Water; DEPC (Diethyl Pyrocarbonate treated Sterilized distilled Water).


## PREREQUISITE

- Applied Biosystems 7500 Fast (ABF), CFX96 (CFX), QuantStudio 5 (QS5), qTOWER<sup>3</sup> (QT3) and Gentier 96E (GTE) Real-time PCR Instrument system.
- Pipettes and Pipettes tips with aerosol barrier.
- Vortex mixer.
- Centrifuge.
- Disposable powder-free gloves.

## WARNING AND PRECAUTION

- This product is designed *In vitro diagnostics* (IVD).
- The test has to be performed by well-trained and qualified personnel.
- Dispose of unused reagents and waste in accordance with country, federal, state and local regulations.
- Unnecessary repeated freezing and thawing will be occurred inaccurate results.
- Do not mix reagent from different batches of the kit.
- Do not modify the reagent/sample volume used in the test or use in a wrong way which is not recommended.

## WORK FLOW



## PROTOCOL

### A. Specimen

This product must be used with viral RNA samples extracted from human respiratory specimens such as Nasopharyngeal & Oropharyngeal (Throat) swab, Sputum and Bronchoalveolar Lavage Fluid (BAL).

### B. RNA Extraction

It is recommended to use commercialized extraction kit such as QIAamp Viral RNA mini kit (Qiagen).

### C. Reagent Preparation

Before setting up PCR, all components need to be thawed, gently mixed and centrifuged briefly to collect solution at the bottom.

- Mix 4 µL of DEPC Water, 4 µL of COVID-19 Fast Reaction Buffer, 5 µL of COVID-19 Fast Probe Mixture and 2 µL of COVID-19 Fast Enzyme Mixture to prepare master mixture per each reaction (refer to the below). Prepare enough volume of master mixture for all the reactions Fast extra to prevent pipetting error.
- After mixing well, place 15 µL of master mixture into 96-well plate or PCR tube.
- Add 5 µL of extracted RNA sample into 96-well plate or tube, then mix all components by pipetting. Proceed in the same way with other RNA samples, Positive Control and Negative Control.
- Accurately close the tube with the cap or seal the 96-well plate.
- Transfer the tubes or 96-well plate for test into the Real-time PCR and start for the amplification.

Component	Per reaction (µL)
DEPC Water	4 µL
COVID-19 Fast Reaction Buffer	4 µL
COVID-19 Fast Probe Mixture	5 µL
COVID-19 Fast Enzyme Mixture	2 µL
Extracted RNA or PC <sup>1)</sup> or NC <sup>2)</sup>	5 µL
<b>Total Volume</b>	<b>20 µL</b>

<sup>1)</sup> PC, Positive Control; <sup>2)</sup> NC, Negative Control

## D. Setting of Real-time PCR

- Referring to the instrument manual, set on the dedicated software for the parameters of thermal cycle.
- Set up the PCR program and fluorescenc as following, and then click the start "Run" button.

PCR program					
Cycle	Temp.	Time			
		ABF	QS5	CFX/GTE	QT3 (Ramping rate, °C/S)**
1 cycle	50 °C	5 min	5 min	5 min	5 min (4)
1 cycle	95 °C	20 sec	20 sec	20 sec	20 sec (4)
45 cycles	95 °C	1 sec	1 sec	1 sec	1 sec (4)
	58 °C*	30 sec	20 sec	5 sec	5 sec (2)

\* Select "Collect Data" or "plate Read"

\*\* For qTOWER<sup>3</sup>, set the Ramping rate as follows.

Fluorescence Setting					
Target	ABF	QS5	CFX/GTE	QT3	
				Reporter dye	Gain
RdRp gene	FAM	FAM	FAM	FAM	3
N gene	JOE	VIC	VIC	HEX-3	5
E gene	Texas Red	ROX	Texas Red	Texas Red	5
IC	Cy5	Cy5	Cy5	Quasar 670	5
NO Quencher, Reference Dye None					

## E. Analysis Setting

The values of fluorescence emitted by the specific probes and by the specific internal control probe in amplification reactions should be analyzed by the instrument software.

- Click Analysis mode after completion and choose analysis setting from Amplification Plot.
- Click "Edit Default Settings" to set threshold values as shown below.

Target	Threshold					Baseline	
	ABF	QS5	CFX	GTE	QT3	Begin	End
RdRp gene	100,000	75,000	1,000	1,000	3	3	15
N gene	50,000	50,000	500	1,000	3	3	15
E gene	50,000	50,000	500	1,000	3	3	15
IC	10,000	10,000	100	300	3	3	15

- For qTOWER<sup>3</sup>, set the color compensation setting to Standard 1.

## F. Result Interpretation

#	RdRP	N	E	IC	Assay Result
1	≤40	≤40	≤40	≤35	COVID-19 Positive
2	≤40	≤40	U.D*	≤35	COVID-19 Positive
3	≤40	U.D	≤40	≤35	COVID-19 Positive
4	≤40	U.D	U.D	≤35	Repeat the test (COVID-19 Positive if RdRp ≤40)
5	U.D	≤40	≤40	≤35	COVID-19 Positive
6	U.D	≤40	U.D	≤35	Repeat the test (COVID-19 Positive if N ≤40)
7	U.D	U.D	≤40	≤35	Sarbecovirus
8	U.D	U.D	U.D	≤35	Negative
9	U.D	U.D	U.D	U.D	Invalid (re-test)

\*Undetermined (or Not applicable (N/A) or No Ct or '-')

## QUALITY CONTROL

- PC and NC Ct range should be as below:

#	RdRp gene	N gene	E gene	IC
PC	≤35	≤35	≤35	≤35
NC	U.D	U.D	U.D	U.D

## PERFORMANCE

Criteria	Result
Analytical Specificity	26 DNA/RNA samples were tested on the GeneFinder™ COVID-19 Fast RealAmp Kit in order to evaluate the possibility of cross-reactivity. 26 DNA/RNA samples which have no concern with the detection target of the kit were negative. * Cross reactivity : 100% Specificity
Analytical Sensitivity	Serial dilutions (5,2,1,0.2copies/ul) of COVID-19 RNA (3 batches, 20 times of repeat test each) were tested. * Analytical Sensitivity : 1) RdRp gene: 5 copies/ul 2) N gene: 2 copies/ul 3) E gene: 5 copies/ul

Reproducibility	Reproducibility was confirmed with identical standard substances at different conditions; different LOT, place, time and testers by 3 batch testing. Criteria of Reproducibility was CV < 5% of Ct Value.
Freeze / Thaw Safety	Freeze/Thaw safety of GeneFinder™ COVID-19 Fast RealAmp Kit was confirmed by 12 times of Freeze/Thaw repeat test. Criteria of safety was CV < 5% of Ct Value.

## TROUBLE SHOOTING

Problem	Possible Cause	Recommendation
No signal in all samples including positive control	Error in Master mixture preparation	Check the dispensing volume during preparation of master mixture
	Inhibitors added	Repeat the extraction step with new sample
	Probe degradation	Use a new probe reagent
	Positive Control degradation	Use a new positive control
	Omitted components	Verify each component, repeat the RT-PCR mixture preparation
Diverse intensity of fluorescent signals	Instrument setting error	Check the position setting for the positive control on the instruments. Check the Thermal cycle settings on sample instrument
	Pipetting error	Make sure that the equal volume of reactants is added in each tube or plate
	Contamination in the outer surface of PCR tubes or Plate	Wear gloves during the experiment
Weak or no fluorescent signal in samples only	Poor RNA quality	Use recommended kit for RNA extraction and store extracted RNA at -70°C
	Insufficient volume of RNA	Repeat PCR reaction with correct volume of RNA